



Recombinant Arabidopsis thaliana Zinc finger protein CONSTANS-LIKE 7 (COL7)

Product Code	CSB-YP863238DOA
Abbreviation	COL7
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9C9A9
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	≥85% (SDS-PAGE)
Sequence	MVVDVESRTA SVTGEKMAAR GCDACMKRSR ASWYCPADDA FLCQSCDASI HSANHLAKRH ERVRLQSSSP TETADKTTSV WYEGFRRKAR TPRSKSCAFE KLLQIESNDP LVPELGGDED DGFFSFSSVE ETEESLNCCV PVFDPFSDML IDDINGFCLV PDEVNNTTTN GELGEVEKAI MDDEGFMGFV PLDMDLEDLT MDVESLLEEE QLCLGFKEPN DVGVIKEENK VGFEINCKDL KRVKDEDEEEE EEAKCENGGS KDS DREASND KDRKTSFLR LDYGAVISAW DNHGSPWKTG IKPECM LGGN TCLPHVGGY EKLMSSDGSV TRQQGRDGGG SDGEREARVL RYKEKRTRL FSKKIRYEVK KLNAEQRPRI KGRFVKRTSL LT
Source	Yeast
Target Names	COL7
Protein Names	Recommended name: Zinc finger protein CONSTANS-LIKE 7
Expression Region	1-392
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.