



Recombinant Arabidopsis thaliana F-box protein At3g08750 (At3g08750)

Product Code	CSB-EP863245DOA
Abbreviation	At3g08750
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9C9Y4
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	≥85% (SDS-PAGE)
Sequence	MASSKCLLLP SLPFELIEEI LYKIPAESLI RFKSTCKKWY NLITEKRFRMY NHLDHYSRPER FIRTYDQQII DPVTEILSDA LIPDEFRLDY PIYSMVHCDG LMLCTCRKWD NSLAVWNPVL REIKWIKPSV CYLHTDYVGI GYDDNVSRDN YKILKLLGRL PKDDSDPNC EIYEFKSDSW KTLVAKFDWD IDIRCNNGVS VKGKMYWIAK KKEDFTIIRF DFSTETFKEI CVCPTYTLVTR LGCFDGDRLS LLLQGEESQG IEVWMTNKLK DKVVSFSQYF NVTTPDLPAL HLQSCMACPG YSIGKRRNIR VWCEGSVDVD DKSYVIITFY EIDEGGVIKQ IETASYGQFE YDDPFICCYV YVPSLIPIQ
Source	E.coli
Target Names	At3g08750
Protein Names	Recommended name: F-box protein At3g08750
Expression Region	1-369
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.