



# Recombinant Mouse Peroxisomal leader peptide-processing protease (Tysnd1)

<b>Product Code</b>	CSB-EP863575MO
<b>Abbreviation</b>	Tysnd1
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9DBA6
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Mus musculus (Mouse)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MGRQWGPSMR VAEQAGCVVS ASRAGQPDAG SWSCSGVILS RNPGLVLCHG GIFTPLRTG SAALTQTGTA FLPGDSCSDD LRLHVQWGPT AASPAGRADQ ELPGLCTPQC ASLGLPEGAP SRARARPLQP PRPAQLLLLL SCPAFRSHFA RLFGADAVDQ WHFVSSAPDD AVSEEEEEEDQ LRALGWFWALL RVQRGAAAE RRGPVVTVAP LGAVVKGAPL LACGSPFGAF CPDIFLNTLS RGVLSNAAGP LLLTDARCLP GTEGGGVFAA RPAGALVALV AAPLCWKARE WVGLTLLCAA APLLQVARWA LARLHPGSAS LSVLLPPPVDV STPRGLPLRD LGPPWAAAAY LVECGTVWGS GVVVAPRLVV TCRHVAPREA ARVLVHSATP KNVAIWGQVV FATQETSPYD IAVVSLEEEL NGVPTPVPAG HFHEGEPVSV VGFGVFGQAC GPSVTSGILS AVVRVDGSPV MLQTTCAVHG GSSGGPLFSS GSGDLLGIVA SNTRDNNTGA TYPHLNFSIP ITVLQPALKQ YSQTGDLGGL RELDHTTEPV RVVWRLQRPL SEVPRSKL
<b>Source</b>	E.coli
<b>Target Names</b>	Tysnd1
<b>Protein Names</b>	Recommended name: Peroxisomal leader peptide-processing protease EC=3.4.21.- Alternative name(s): Trypsin domain-containing protein 1 Cleaved into the following 2 chains: 1. Peroxisomal leader peptide-processing protease, 10 kDa form
<b>Expression Region</b>	1-568
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final



concentration of glycerol is 50%. Customers could use it as reference.

### Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.