



# Recombinant Rat MARCKS-related protein (Marcksl1)

|                          |   |
|--------------------------|---|
| <b>Product Code</b>      | CSB-MP863634RA  |
| <b>Abbreviation</b>      | Marcksl1  |
| <b>Storage</b>           | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.  |
| <b>Uniprot No.</b>       | Q9EPH2  |
| <b>Product Type</b>      | Recombinant Protein   |
| <b>Immunogen Species</b> | Rattus norvegicus (Rat)   |
| <b>Purity</b>            | >85% (SDS-PAGE)   |
| <b>Sequence</b>          | GSQSSKAPR GDVTAEAAAG ASPAKANGQE NGHVKSNGLD TPKGEGESPP<br>VNGADEAAAGA TGDAIEPAPP SQEAEAKGEV APKETPKKKK KFSFKKPFKL<br>SGLSFKRNRK EGGGDSSASS PTEEEQEQGE ISACSDEGTA QEGKAAATPE<br>SQEPQAKGAE ASAVSKGGDA EEEAGPQAAE PSTPSGPESG PASASAENE   |
| <b>Source</b>            | Mammalian cell  |
| <b>Target Names</b>      | Marcksl1  |
| <b>Protein Names</b>     | Recommended name: MARCKS-related protein Alternative name(s): Brain protein F52 MARCKS-like protein 1 Macrophage myristoylated alanine-rich C kinase substrate Short name= Mac-MARCKS Short name= MacMARCKS   |
| <b>Expression Region</b> | 2-199   |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.   |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.   |
| <b>Protein Length</b>    | Full Length of Mature Protein   |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| <b>Shelf Life</b>        | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.  |