



Recombinant Mouse Vacuolar protein sorting-associated protein 35 (Vps35)

Product Code	CSB-MP863646MO
Abbreviation	Vps35
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9EQH3
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	≥85% (SDS-PAGE)
Sequence	MPTTQQSPQD EQEKLLDEAI QAVKVQSFQM KRCLDKNKLM DALKHASNML GELRTSMLSP KSYEELYMAI SDELHYLEVY LTDEFKGRK VADLYELVQY AGNIIPRLYL LITVGVVYVK SFPQSRKDIL KDLVEMCRGV QHPLRGLFLR NYLLQCTRNI LPDEGEPTDE ETTGDISDSM DFVLLNFAEM NKLWVRMQHQ GHSRDREKRE RERQELRILV GTNLVRLSQL EGVNVERYKQ IVLTGILEQV VNCRDALAEQ YLMECIIQVF PDEFHLQTLN PFLRACAEHL QNVNVKNIII ALIDRLALFA HREDGPGIPA EIKLFDIFSQ QVATVIQSRQ DMPSEDVVS LQVSLINLAMK CYPDRVDYVD KVLETTVEIF NKLNLEHIAT SSAVSKELTR LLKIPVDYTN NILTVLKLKH FHPLFEYFDY ESRKSMSCYV LSNVLDYNT EIVSQDQVDSI MNLVSTLIQD QPDQPVDPD PEDFADEQSL VGRFIHLLRS DDPDQQYLIL NTARKHFGAG GNQRIRFTLP PLVFAAYQLA FRYKENSQMD DKWEKQCQKI FSFAHQTISA LIKAEALP LRLFLQGALA AGEIGFENHE TVAYEFMSQA FSLYEDEISD SKAQLAAITL IIGTFERMKC FSEENHEPLR TQCALAASKL LKKPDQGRAV STCAHLFWSG RNTDKNGEEL HGGKRVMECL KKALKIANQC MDPSLQVQLF IEILNRYIYF YEKENDAVTI QVLNQLIQKI REDLPNLESS EETEIQINKHF HNTLEHLRSR RESPESEGPI YEGLIL
Source	Mammalian cell
Target Names	Vps35
Protein Names	Recommended name: Vacuolar protein sorting-associated protein 35 Alternative name(s): Maternal-embryonic 3 Vesicle protein sorting 35
Expression Region	1-796
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Target Details	This gene belongs to a group of vacuolar protein sorting (VPS) genes. The encoded protein is a component of a large multimeric complex, termed the



retromer complex, involved in retrograde transport of proteins from endosomes to the trans-Golgi network. The close structural similarity between the yeast and human proteins that make up this complex suggests a similarity in function. Expression studies in yeast and mammalian cells indicate that this protein interacts directly with VPS35, which serves as the core of the retromer complex.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.