



Recombinant Rat Platelet-derived growth factor C (Pdgfc)

Product Code	CSB-EP863652RA-B
Abbreviation	Pdgfc
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9EQX6
Product Type	Recombinant Protein
Immunogen Species	Rattus norvegicus (Rat)
Purity	>85% (SDS-PAGE)
Sequence	ESNLSSKL QLSSDKEQNG VQDPRHERVV TISGNNGSIHS PKFPHTYPRN TVLVWRLVAV DENVRIQLTF DERFGLEDPE DDLCKYDFVE VEPPSDGSLV GRWCGSGTVP GKQTSKGNHI RIRFVSDEYF PSEPGFCIH YSIIMPQVTET TSPSVLPSSA LSLDLLNNAV TAFSTVEELI RFLEPDRWQI DLDSLYKPTW PLLGKAFLYG KKS KAVNLNL LKEEVKLYSC TPRNFSVSIR EELKRTDTIF WPGCLLVKRC GGNCACCLHN CNECQCVP RK VTKKYHEVLQ LRPKIGVKGL HKSLTDVALE HHEECDCVCR GNTEG
Source	E.coli
Target Names	Pdgfc
Protein Names	Recommended name: Platelet-derived growth factor C Short name= PDGF-C Alternative name(s): Fallotein Spinal cord-derived growth factor Short name= rScdfg VEGF-E Cleaved into the following 2 chains: 1. Platelet-derived gro
Expression Region	23-345
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Target Details	This protein is a member of the platelet-derived growth factor family. The four members of this family are mitogenic factors for cells of mesenchymal origin and are characterized by a core motif of eight cysteines. This gene product appears to form only homodimers. It differs from the platelet-derived growth factor alpha and beta polypeptides in having an unusual N-terminal domain, the CUB domain.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final



concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.