



Recombinant Arabidopsis thaliana Transcription factor MYB75 (MYB75)

Product Code	CSB-EP863725DOA
Abbreviation	MYB75
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9FE25
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	≥85% (SDS-PAGE)
Sequence	MEGSSKGLRK GAWTTEEDSL LRQCINKYGE GKWHQVPVRA GLNRCRKSCR LRWLNLYLKPS IKRGKLSSDE VDLLLLRLHRL LGNRWSLIAG RLPGRTANDV KNYWNTHL SK KHEPCCKIKM KKRDIPIPT TPALKNNVYK PRPRSFTVNN DCNHLNAPPK VDVNPPCLGL NINNVCDNSI IYNKDKKKDQ LVNNLIDGDN MWLEKFLEES QEVDILVPEA TTTEKGDTLA FDVDQLWSLF DGETVKFD
Source	E.coli
Target Names	MYB75
Protein Names	Recommended name: Transcription factor MYB75 Alternative name(s): Myb-related protein 75 Short name= AtMYB75 Production of anthocyanin pigment 1 protein Suc-induced anthocyanin accumulation locus 1
Expression Region	1-248
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.