



Recombinant Arabidopsis thaliana Probable UDP-arabinopyranose mutase 5 (RGP5)

Product Code	CSB-YP863733DOA
Abbreviation	RGP5
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9FFD2
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	≥85% (SDS-PAGE)
Sequence	MSLAIEINKNE VDIVIGALNA DLTQFLTSWR PFFSGFHLIV VKDPELKEEL NIEGFDVDV YSKTDMKVV GASNSTMFSG YSCRYFGYLV SKKKYIVSID DDCVPAKDPK GFLVDAVTQH VINLENPATP LFFNTLYDPY CEGADFVRGY PFSLRSGVPC AASCGLWLNL ADLDAPTQAL KTEKRNTAYV DAVMTVPAKA MLPISGINIA FNRELVGPAL VPALRLAGEG KVRWETLEDV WCGMCLKHIS DHLGYGVKTG LPYVWRNERG DAVESLRKKW EGMKLMEKSV PFFDSLKLPE TALKVEDCVI ELAKAVKEQL GSDDPAFTQA ADAMVKWVQL WNSVNSSA
Source	Yeast
Target Names	RGP5
Protein Names	Recommended name: Probable UDP-arabinopyranose mutase 5 EC= 5.4.99.30 Alternative name(s): Reversibly glycosylated polypeptide 5 Short name= AtRGP5 UDP-L-arabinose mutase 5
Expression Region	1-348
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.