



Recombinant Arabidopsis thaliana FBD-associated F-box protein At5g38590 (At5g38590)

Product Code	CSB-MP863738DOA
Abbreviation	At5g38590
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9FFW2
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	≥85% (SDS-PAGE)
Sequence	MDKINGLPDD LLVKILSYVP TDIKVSTSL SKRWEFLWMW LPNLDYTSRW CRKPGDVGLR DFIHKNLPLH RAPVIESLRF HSNSPDIKPE DIRRWIEIAV SRHVHDLID HFSEENIFL SSFFACKSLV TLKLRVTLR DIPSMVCLPS LKTTTTDNVS FVEGKSLQEL LSICPVLEDL SVYCDDYENT KELITVPSL LSLSLYIPDE WLLDGYWIDT PSLEYLKLED WNSCDHLSLI KNMPKLREAY VDAKCFLPKS VIESITSVKH LTICSKDGYG DGFVFNQLEH LTLVCVRGDS PSLLGQLLKD SPNLRILEIS VMEDHVDDVG ISLDGWNQPS SVPECLLSSL QIFKWPQYLG RPEDRDIAVY ILKNARHLKK TTILADRCER FVTQRRMIKE LLQALPARIC
Source	Mammalian cell
Target Names	At5g38590
Protein Names	Recommended name: FBD-associated F-box protein At5g38590
Expression Region	1-410
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.