



Recombinant Arabidopsis thaliana F-box protein At5g49610 (At5g49610)

Product Code	CSB-EP863750DOA-B
Abbreviation	At5g49610
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9FGY4
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MDNQKGALFP DEVILQILAR LPVKSLFRFK SVCKSWYRLP SDKYFTSLFN QLSVKEQLLV AEVSDSSSLI CVDNLRGVSE LSLDFVRDRV RIRVSSNGLL CCSSIPEKGV YYVCNPSTRE YRKLPKSRER PVTRFYDPGE ATLVGLACDL SKNKFNVVLA GYHRSFGQRP DGSFICLVFD SESNKWRKFV SVLEECSTFH MSKNQVVFVN GMLHWLMSG L CYILALDVEH DVWRKISLPD EIRIGNGGGN RVYLLESDGF LSVIQLSDVW MKIWKMSEYE TETWSVVD SI SLRCIKGLVP GIFPICQTGE YVFLATHKQV L VYQRRSKLW KEMFSVKGSS SLPLWFSAHA FRSTIVPCN
Source	E.coli
Target Names	At5g49610
Protein Names	Recommended name: F-box protein At5g49610
Expression Region	1-359
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.