



Recombinant Arabidopsis thaliana Protein ABIL4 (ABIL4)

Product Code	CSB-MP863759DOA
Abbreviation	ABIL4
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9FHY1
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	≥85% (SDS-PAGE)
Sequence	MASSTSLAIV LHQSSNHDEL FMKQTLQFSE TLKDLKnlrk QLYSAAEYFE TSYGKAEHKE TVIETLKEYA AKAVVNTVDH LGSVSDKFNS FLSDNSTHFS TTHLRLSSLE QRMRLCRDYM GKSGTHQHLL LFQYPRHHKR YFFPQQGRGT SFSAGDDSHR FTS AVRSTIL ENLPNTARKA NKTGSF SFAP IVHNNINNRT PNKRSNSPMR FPLLRSGSL KRSSSPSQPK KPPLALPEPQ RAISVSRNTE IVEIKQSSSR KGKKILMFKA LMSMSKSRN
Source	Mammalian cell
Target Names	ABIL4
Protein Names	Recommended name: Protein ABIL4 Alternative name(s): Abl interactor-like protein 4 Short name= AtABIL4
Expression Region	1-279
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.