



Recombinant Arabidopsis thaliana F-box/FBD/LRR-repeat protein At5g56420 (At5g56420)

Product Code	CSB-EP863810DOA-B
Abbreviation	At5g56420
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9FM89
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	≥85% (SDS-PAGE)
Sequence	MVSYRDRLSQ LPDDFLLQIL SWLPTKDVLV TSLLSKRWRF LWTLVPRLNY DLRLHDNTCP RFSQFVDRSL LLHKAPTLES LNIKIGSICF TAEKDVGWVW RIGVDRFVRE LSVSYCSGEE PIRLPKCLFT CSTLAVLKLE NITLEDASCY VCFQSLKTLH LLDVKYLLDDQ SLPRISSCS SLEDLVVQRC PGDNVKVVTV TAPSLKTLSL HKSSQAFEGD DDGFLIDTPK LKRVDIEDYW GGFCYIENMP EVVEANVDVI YKNTEKLLGS ITSVKRLALC LITSDAAYPA GTIFSQVLHL ELCTCAPRWW DLLTRLIEDS PKLRVCLKRQ KHIRRAPSPR ASWKQPALPK CLLFHLETFK WELYEGSQKQ KEVATFILKH AIRLKTAIIS PKPTSTLLEK HEMLKDLSSS PRGSSTCELL FD
Source	E.coli
Target Names	At5g56420
Protein Names	Recommended name: F-box/FBD/LRR-repeat protein At5g56420
Expression Region	1-422
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.