



Recombinant Arabidopsis thaliana FBD-associated F-box protein At5g22730 (At5g22730)

Product Code	CSB-EP863825DOA-B
Abbreviation	At5g22730
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9FNI8
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	≥85% (SDS-PAGE)
Sequence	<p> MLYEGYLKYM KAKYGERSQG TYFMAGEDLI SKLPDSLITQ ILLYLPIKDI VRTSSLSSRW KSLWLLIPRL DLDSEEFQDY NAFVGFMNKF IDFSGEEKIC LDKLKLSSRK TVNDLPCVTR WIDFVVRK KLHLDVECLVN RKFLEEMPLS LYVCDTLVNL RLHRVLLGKF EAVSLPCLKT MRLEENVYAN DVVLESLISS CPVLKDLIIL RMFEDNVKVL RVHSLTLTSL NIDFNFGEGD DFVDGFDKKV SGVLIDAPRL KYLKFQDDL SSKIITNSGS LAKVNVVYVF NENDCADVVD IPRRNMVRNF LTSISGVSDM KISQHFVEFL YYYKDFDPLP QFCNLSRLKA EISLYFLEIL PTILESCP NLKSLVMVLEFY LQEEDEPIIF SSVPRCLVSS LESVEIKKFN GRPAKMEVAR YFLENSGVLQ KLVHLRCST HEDGFYILKD LLALPRGSST CRIVVC </p>
Source	E.coli
Target Names	At5g22730
Protein Names	Recommended name: FBD-associated F-box protein At5g22730
Expression Region	1-466
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.