



# Recombinant Arabidopsis thaliana FBD-associated F-box protein At5g22730 (At5g22730)

<b>Product Code</b>	CSB-EP863825DOA
<b>Abbreviation</b>	At5g22730
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9FNI8
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Arabidopsis thaliana (Mouse-ear cress)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	<p>           MLYEGYLKYM KAKYGERSQG TYFMAGEDLI SKLPDSLITQ ILLYLPIKDI            VRTSSLSSRW KSLWLLIPRL DLDSEEFQDY NAFVGFMNKF IDFSGEEKIC            LDKLKLSSRK TVNDLPCVTR WIDFVVRK KLHLDVECLVN RKFLEEMPLS            LYVCDTLVNL RLHRVLLGKF EAVSLPCLKT MRLEENVYAN DVVLESLISS            CPVLKDLIIL RMFEDNVKVL RVHSLTLTSL NIDFNFGEGD DFVDGFDKKV            SGVLIDAPRL KYLKFQDDLS GSKIITNSGS LAKVNVVYVF NENDCADVVD            IPRRNMVRNF LTSISGVSDM KISQHFVEFL YYYKDFDPLP QFCNLSRLKA            EISLYFLEIL PTILESCP NL KSLVMVLEFY LQEEDEPIIF SSVPRCLVSS            LESVEIKKFN GRPAKMEVAR YFLENSGVLQ KLVLHLRCST HEDGFYILKD            LLALPRGSST CRIVVC         </p>
<b>Source</b>	E.coli
<b>Target Names</b>	At5g22730
<b>Protein Names</b>	Recommended name: FBD-associated F-box protein At5g22730
<b>Expression Region</b>	1-466
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.