



Recombinant Arabidopsis thaliana Agamous-like MADS-box protein AGL31 (AGL31)

Product Code	CSB-MP863830DOA
Abbreviation	AGL31
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9FPN7
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MGRKKVEIKR IENKSSRQVT FSKRRNGLIE KARQLSILCE SSI AVL VVSG SGKLYKSASG DNMSKIIDRY EIHHAELEA LD LAEKTRNY LPLKELLEIV QSKLEESNVD NASVDTLISL EEQLETALSV TRARKTE LMM GEV KSLQKTE NLLREENQTL ASQV GKK TFL VIEGDRGMSW ENGSGNKVRE TLPLLK
Source	Mammalian cell
Target Names	AGL31
Protein Names	Recommended name: Agamous-like MADS-box protein AGL31 Alternative name(s): Protein MADS AFFECTING FLOWERING 2
Expression Region	1-196
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.