



Recombinant Arabidopsis thaliana Putative F-box protein At1g67623 (At1g67623)

Product Code	CSB-EP863856DOA-B
Abbreviation	At1g67623
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9FXC7
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	≥85% (SDS-PAGE)
Sequence	MVVSLSLHAS SSIRKSKTKA SLCLDSLPEL LLVEISSCTG ASSLSAVRNL RLVSKSFRRRI CDEKYVYFRL SLKEIEFLPW HENSAKFIER CTESRNPEAL FQKGFINYFR DKLQDRGLEY LAEAAEKGIK EAKYVYGVIL ICLGGKTKQK GFEILSSVIK QLMSTTMNEL VEFYKIQKI RYGFWWSNT VVEQLKTAYV SEKCKDCKT RMLLLVMNRG WYLFGEETDL DFSSACELYL YINRNFPEA KKS KIDGEIY WACEFEPNRF VRTEPNRGLD KSVRLRLIEY NRSVSV
Source	E.coli
Target Names	At1g67623
Protein Names	Recommended name: Putative F-box protein At1g67623
Expression Region	1-296
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.