



Recombinant Inosine triphosphate pyrophosphatase (hap-1)

Product Code	CSB-YP863924CXY
Abbreviation	hap-1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9GYG4
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	≥85% (SDS-PAGE)
Sequence	MSLRKINFVT GNVKKLEEVK AILKNFEVSN VDVLDDEFQG EPEFIAERKC REAVEAVKGP VLVEDTSLCF NAMGGLPGPY IKWFLKNLKP EGLHNMLAGF SDKTAYAQCI FAYTEGLGKP IHVFAGKCPG QIVAPRGDTA FGWDPCFQPD GFKETFGEMD KDVKNEISHR AKALELLKEY FQNN
Source	Yeast
Target Names	hap-1
Protein Names	Recommended name: Inosine triphosphate pyrophosphatase Short name= ITPase Short name= Inosine triphosphatase EC= 3.6.1.19 Alternative name(s): Non-canonical purine NTP pyrophosphatase Non-standard purine NTP pyrophosphatas
Expression Region	1-184
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.