



Recombinant Human PR domain zinc finger protein 14 (PRDM14)

Product Code	CSB-EP863934HU-B
Abbreviation	PRDM14
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9GZV8
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MALPRPSEAV PQDKVCYPPE SSPQNLAAYY TPFPSYGHYR NSLATVEEDF QPFRQLEAAA SAAPAMPPFP FRMAPPLLSP GLGLQREPLY DLPWYSKLPP WYPIPHVPRE VPPFLSSSHE YAGASSEDLG HQIIGGDNES GPCCGPTLI PPPPADASLL PEGLRTSQLL PCSPSKQSED GPKPSNQEKG SPARFQFTEE DLHFVLYGVT PSLEHPASLH HAISGLLVPP DSSGSDSLPQ TLDKDSLQLP EGLCLMQTVF GEVPHFGVFC SSFIKGVRF GPFQKGKVVNA SEVKTYGDNS VMWEIFEDGH LSHFIDGKGG TGNWMSYVNC ARFPKEQNLV AVQCQGHIFY ESCKEIQNQ ELLVWYGDCY EKFLDIPVSL QVTEPGKQPS GPSEESAEGY RCERCGKVFT YKYIRDKHLK YTPCVDKGDR KFPCSLCKRS FEKRDLRIH ILHVHEKHRP HKCSTCGKCF SQSSSLNKHM RVHSGDRPYQ CVYCTKRFTA SSILRTHIRQ HSGEKPFKCK YCGKSFASHA AHDShVRRSH KEDDGCSCSI CGKIFSDQET FYSHMKFHED Y
Source	E.coli
Target Names	PRDM14
Protein Names	Recommended name: PR domain zinc finger protein 14 Alternative name(s): PR domain-containing protein 14
Expression Region	1-571
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.