



# Recombinant Human Cysteine/serine-rich nuclear protein 2 (CSRNP2)

<b>Product Code</b>	CSB-BP863949HU
<b>Abbreviation</b>	CSRNP2
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9H175
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MDAFTGSGLK RKFDDVDVGS SVSNSDDEIS SSSDADSCDS LNPPTTASFT PTSILKRQKQ LRRKNVRFQDQ VTVYYFARRQ GFTSVPSQGG SSLGMAQRHN SVRSYTLCEF AQEQEVNHRE ILREHLKEEK LHAKMKLTK NGTVESVEAD GLTLDDVSDE DIDVENVEVD DYFFLQPLPT KRRRALLRAS GVHRIDAEK QELRAIRLSR EECGCDCRLY CDPEACACSQ AGIKCQVDRM SFPCGCSR DG CGNMAGRIEF NPIRVRTHYL HTIMKLELES KRQVSRPAAP DEEPSPTASC SLTGAQGSET QDFQEFIAEN ETAVMHLQSA EELERLKAEE DSSGSSASLD SSIESLGVCI LEEPLAVPEE LCPGLTAPIL IQAQLPPGSS VLCFTENSDH PTASTVNSPS YLNSGPLVYY QVEQRPVLGV KGEPGTEEGS ASFPKEKDLN VFSLPVTSLV ACSSTDPAAL CKSEVGKTPPT LEALLPEDCN PEEPENEDFH PSWSPSSLPF RTDNEEGCGM VKTSQQNEDR PPEDSSLELP LAV
<b>Source</b>	Baculovirus
<b>Target Names</b>	CSRNP2
<b>Protein Names</b>	Recommended name: Cysteine/serine-rich nuclear protein 2 Short name= CSRNP-2 Alternative name(s): Protein FAM130A1 TGF-beta-induced apoptosis protein 12 Short name= TAIP-12
<b>Expression Region</b>	1-543
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.