



Recombinant Human Mediator of RNA polymerase II transcription subunit 28 (MED28)

Product Code	CSB-EP863955HU
Abbreviation	MED28
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9H204
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MAAPLGGMFS GQPPGPPQAP PGLPGQASLL QAAPGAPRPS SSTLVDELES SFEACFASLV SQDYVNGTDQ EEIRTGVDQC IQKFLDIARQ TECFFLQKRL QLSVQKPEQV IKEDVSELRN ELQRKDALVQ KHLTKLRHWQ QVLEDINVQH KKPADIPQGS LAYLEQASAN IPAPLKPT
Source	E.coli
Target Names	MED28
Protein Names	Recommended name: Mediator of RNA polymerase II transcription subunit 28 Alternative name(s): Endothelial-derived protein 1 Mediator complex subunit 28 Merlin and Grb2-interacting cytoskeletal protein Short name= Magicin Tumor a
Expression Region	1-178
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.