



Recombinant Human Cytosolic Fe-S cluster assembly factor NARFL (NARFL)

Product Code	CSB-EP863991HU
Abbreviation	NARFL
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9H6Q4
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	ASPFSGALQ LTDLDDFIGP SQECIKPVKV EKRAGSGVAK IRIEDDGSYF QINQDGGTRR LEKAKVSLND CLACSGCITS AETVLITQQS HEELKKVLDA NKMAAPSQQR LVVVSVSPQS RASLAARFQL NPTDTARKLT SFFKKIGVHF VFDTAFSRHF SLLESQREFV RFRGQADCR QALPLLASAC PGWICYAEKT HGSFILPHIS TARSPQQVMG SLVKDFFAQQ QHLTPDKIYH VTVMPCYDKK LEASRPDFFN QEHQTRDVC VLTGGEVFRLL LEEEGVSLPD LEPAPLDSL SGASAEPTS HRGGGSGGYL EHVFRHAARE LFGIHVAEVT YKPLRNKDFQ EVTLEKEGQV LLHFAMAYGF RNIQNLVQRL KRGRCPYHYV EVMACPSGCL NGGGQLQAPD RPSRELLQHV ERLYGMVRAE APEDAPGVQE LYTHWLQGTD SECAGRLLHT QYHAVEKAST GLGIRW
Source	E.coli
Target Names	NARFL
Protein Names	Recommended name: Cytosolic Fe-S cluster assembly factor NARFL Alternative name(s): Iron-only hydrogenase-like protein 1 Short name= IOP1 Nuclear prelamins A recognition factor-like protein Protein related to Narf
Expression Region	2-476
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.