



Recombinant Human 5-azacytidine-induced protein 2 (AZI2)

Product Code	CSB-MP863992HU
Abbreviation	AZI2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9H6S1
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MDALVEDDIC ILNHEKAHKR DVTVPVSIYS GDESVASHFA LVTAYEDIKK RLKDSEKENS LLKKRIRFLE EKLIARFEEE TSSVGREQVN KAYHAYREVC IDRDNLKSKL DKMNKDENSES LKVLNEQLQS KEVELLQLRT EVETQQVMRN LNPPSSNWEV EKLSCDLKIH GLEQELELMR KECSDLKIEL QKAKQTDPYQ EDNLKSRDLQ KLSISSDNMQ HAYWELKREM SNLHLVTQVQ AELLRKLKTS TAIKKACAPV GCSEDLGRDS TKLHLMNFTA TYTRHPPLLP NGKALCHTTS SPLPGDVKVL SEKAILQSWT DNERSIPNDG TCFQEHSSYG RNSLEDNSWV FPSPPKSSET AFGETKTKTL PLPNLPPLHY LDQHNQNCLY KN
Source	Mammalian cell
Target Names	AZI2
Protein Names	Recommended name: 5-azacytidine-induced protein 2 Alternative name(s): NF-kappa-B-activating kinase-associated protein 1 Short name= Nak-associated protein 1 TILP
Expression Region	1-392
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.