



# Recombinant Human Tudor domain-containing protein 3 (TDRD3)

<b>Product Code</b>	CSB-YP863996HU
<b>Abbreviation</b>	TDRD3
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9H7E2
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MLRLQMTDGH ISCTAVEFSY MSKISLNTTP GTKVKLSGIV DIKNGFLLL DSNTTVLGGE VEHLIEKWEL QRSLSKHNRS NIGTEGGPPP FVPFGQKCVS HVQVDSRELD RRKTLQVTMP VKPTNDNDEF EKQRTAAIAE VAKSKETKTF GGGGGGARSN LNMNAAGNRN REVLQKEKST KSEGKHEGVY RELVDEKALK HITEMGFSKE ASRQALMDNG NNLEAALNVL LTSNKQKPVM GPPLRGRGKG RGRIRSEDEE DLGNARPSAP STLFDFLESK MGTLNVEEPK SQPQQLHQGQ YRSSNTEQNG VKDNNHLRHP PRNDTRQPRN EKPPRFQRDS QNSKSVLEGS GLPRNRGSER PSTSSVSEVW AEDRIKCDRP YSRYDRTKDT SYPLGSQHSD GAFKKRDNSM QSRSGKGPSF AEAKENPLPQ GSVDYNNQKR GKRESQTSIP DYFYDRKSQT INNEAFSGIK IEKHFVNNTD YQNPVRSNSF IGVPNGEVEM PLKGRRIGPI KPAGPVTAVP CDDKIFYNSG PKRRSGPIKP EKILESSIPM EYAKMWKPGD ECFALYWEDN KFYRAEVEAL HSSGMTAVVK FIDYGNYYEV LLSNIKPIQT EAWEEEGTYD QTLEFRRGGD GQPRRSTRPT QQFYQPPRAR N
<b>Source</b>	Yeast
<b>Target Names</b>	TDRD3
<b>Protein Names</b>	Recommended name: Tudor domain-containing protein 3
<b>Expression Region</b>	1-651
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.