



Recombinant Human Phosphorylated adapter RNA export protein (PHAX)

Product Code	CSB-YP863999HU
Abbreviation	PHAX
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9H814
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	ALEVGD MED GQLSDSDSDM TVAPSDRPLQ LPKVLGGDSA MRAFQNTATA CAPVSHYRAV ESVDSSSEESF SDSDDSDCLW KRKRQKCFNP PPKPEPFQFG QSSQKPPVAG GKKINNIWGA VLQEQNQDAV ATELGILGME GTIDRSRQSE TYNYLLAKKL RKESQEHTKD LDKELDEYMH GGKKMGSKKEE ENGQGHLLKRK RPVKDRLG NR PEMNYKGRYE ITAEDSQEKV ADEISFRLQE PKKDLIARVV RIIGNKKAIE LLMETAEEVEQ NGGLFIMNGS RRRTPGGVFL NLLKNTPSIS EEQIKDIFYI ENQKEYENKK AARKRRTQVL GKKMKQAIKS LNFQEDDDTS RETFASDTNE ALASLDESQE GHAEAKLEAE EAIEVDHSHD LDIF
Source	Yeast
Target Names	PHAX
Protein Names	Recommended name: Phosphorylated adapter RNA export protein Alternative name(s): RNA U small nuclear RNA export adapter protein
Expression Region	2-394
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.