



Recombinant Human WD repeat-containing protein 41 (WDR41)

Product Code	CSB-MP864011HU
Abbreviation	WDR41
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9HAD4
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MLRWLIGGGR EPQGLAEKSP LQTIGEEQTQ NPYTELLVLK AHHDIVRFLV QLDDYRFASA GDDGIVVWN AQTGEKLEL NGHTQKITAI ITFPSLESCE EKNQLITAS ADRTVIVWDG DTTRQVQRIS CFQSTVKCLT VLQRLDVWLS GGNDLCVWNR KLDLLCKTSH LSDTGISALV EIPKNCVVAA VGKELIIFRL VAPTEGSLEW DILEVKRLLD HQDNILSLIN VNDLSFVTGS HVGELIIWDA LDWTMQAYER NFWDPSPQLD TQQEIKLCQK SNDISIHFT CDEENVFAAV GRGLYVYSLQ MKRVIACQKT AHDSNVLHVA RLPNRQLISC SEDGSVRIWE LREKQQLAAE PVPTGFFNMW GFGRVSKQAS QPVKKQQENA TSCSLELIGD LIGHSSSVEM FLYFEDHGLV TCSADHLIIL WKNGERESGL RSLRLFQKLE ENGDLYLAV
Source	Mammalian cell
Target Names	WDR41
Protein Names	Recommended name: WD repeat-containing protein 41
Expression Region	1-459
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.