



Recombinant Human RNA polymerase II elongation factor ELL3 (ELL3)

Product Code	CSB-MP864017HU
Abbreviation	ELL3
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9HB65
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MEELQEPLRG QLRLCFTQAA RTSLLLLRLN DAALRALQEC QRQQVRPVIA FQGHRYLRL PPGWSCFLS FIVSQCCQEG AGGSLDLVCQ RFLRSGPNSL HCLGSLRERL IIWAAMDSIP APSSVQGHNL TEDARHPESW QNTGGYSEGD AVSQPQMALE EVSVSDPLAS NQGQSLPGSS REHMAQWEVR SQTHVPNREP VQALPSSASR KRLDKKRSVP VATVELEEKR FRTLPLVPSP LQGLTNQDLQ EGEDWEQEDE DMDPRLEHSS SVQEDSESPS PEDIPDYLLQ YRAIHSAEQQ HAYEQDFETD YAEYRILHAR VGTASQRFIE LGAEIKRVRR GTPEYKVLED KIIQEYKKFR KQYPSYREEK RRCEYLHQKL SHIKGLILEF EEKNRGS
Source	Mammalian cell
Target Names	ELL3
Protein Names	Recommended name: RNA polymerase II elongation factor ELL3
Expression Region	1-397
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.