



Recombinant Human Putative N-acetylated-alpha-linked acidic dipeptidase (FOLH1B)

Product Code	CSB-BP864018HU
Abbreviation	FOLH1B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9HBA9
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MGGSAPPDSS WRGSLKVSYN VGPGFTGNFS TQKVKMHIHS TNEVTRIYNV IGTLRGAVEP DRYVILGGHR DSWVFGGIDP QSGAAVVHET VRSFGTLKKE GWRPRRILF ASWDAEEFGL LGSTEWAEDN SROLLQERGVA YINADSSIEG NYTLRV DCTP LMYSLVYNLT KELKSPDEGF EGKSLYESWT KKSPSPEFSG MPRISKLGSG NDFEVFFQRL GIASGRARYT KNWETNKFSG YPLYHSVYET YELVEKFYDP MFKYHLTVAQ VRGGMVFELA NSIVLPFDCR DYAVVLRKYA DKIYNISMKH PQEMKTYSL S FDSLFSAVKN FTEIASKFSE RLQDFDKSNP ILLRMMNDQL MFLERAFIDP LGLPDRPFYR HVIYAPSSHN KYAGESFPGI YDALFDIESK VDPSKAWGDV KRQISVAAFT VQAAAETLSE VA
Source	Baculovirus
Target Names	FOLH1B
Protein Names	Recommended name: Putative N-acetylated-alpha-linked acidic dipeptidase Short name= NAALADase EC= 3.4.-.- Alternative name(s): Cell growth-inhibiting gene 26 protein Prostate-specific membrane antigen-like protein Putative fo
Expression Region	1-442
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.