



Recombinant Danio rerio Hematopoietically-expressed homeobox protein hhex (hhex)

Product Code	CSB-EP864272DIL
Abbreviation	hhex
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9IAV3
Product Type	Recombinant Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Purity	>85% (SDS-PAGE)
Sequence	MQFQHPHAPL YAPAPAPPAH PTPFYIEDIL GRTGSSSSGPV VPTPTLPSPN SSFTSLIPSY RTPIYELTPI HPVLSQYASM YPFQRSVGDF AHALIRHDPL GKPLLWSPFI QRPLHKRKGQ QVRFSDNDQTI ELEKKFETQK YLSPPERKRL AKMLQLSERQ VKTWFQNRRA KWRRLKQENP PSTGKREAED SDTRRLSDAA ARARELESGA STDSEELLDI EDEHQFTL
Source	E.coli
Target Names	hhex
Protein Names	Recommended name: Hematopoietically-expressed homeobox protein hhex Short name= Homeobox protein hex
Expression Region	1-228
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.