



Recombinant Mouse Coiled-coil domain-containing protein 86 (Ccdc86)

Product Code	CSB-BP864335MO
Abbreviation	Ccdc86
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9JJ89
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	>85% (SDS-PAGE)
Sequence	MDTPLRRSRR LEGLKPLSPE NLPVPEVSRA KRALVDFKSN SEETGELKST RVPLSLPSP GPQPETSPGS PCPPLSLPSP GPQPETSPGS PCPPLSLPSP GPQPETSPGS PCPPLSLPSP GPQPETSPGS PCPPLSLPSP GPQPEASPGS PGPRQDADDG SPQRQPEPHP GSLQPHQDLG LESPAGQTES SPESPQREQP SKLPPPQGEL DSEAAHAKEE VIPGSPEPCP GQQAPGPEPS QPAQELTVQA PSSPERQLEP GKLPPAGETV TESLNLKRV IASPQAPASK KLKEKEELPV IPKGKPKSGR VWKDRSKKRF SQMVQDKPLR TSWQRKMKER QERKLAK DFA RHLEEEKQRR RQEKERRAE NLRRLENER KAEIVQVIRN PAKLKAKKK QLRSEKRDT LALLQKQPPQ RPVAKV
Source	Baculovirus
Target Names	Ccdc86
Protein Names	Recommended name: Coiled-coil domain-containing protein 86 Alternative name(s): Cytokine-induced protein with coiled-coil domain Short name= mCyclon
Expression Region	1-426
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.