



Recombinant Arabidopsis thaliana ATP sulfurylase 1, chloroplastic (APS1)

Product Code	CSB-MP864769DOA
Abbreviation	APS1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9LIK9
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	AGL IAPDGGKLVE LIVEEPKRRE KKHEAADLPR VELTAIDLQW MHVLSSEGWAS PLGGFMRESE FLQTLHFNSL RLDDGSVVM SVPIVLAIDD EQKARIGEST RVALFNSDGN PVAILSDIEI YKHPKEERIA RTWGTTAPGL PYVDEAITNA GNWLIGGDLE VLEPVKYNDG LDRFRLSPAELRKELEKRNA DAVFAFQLRN PVHNGHALLM TDTRRRLLEM GYKNPILLH PLGGFTKADD VPLDWRMKQH EKVLEDGVLD PETTVVSIFP SPMHYAGPTE VQWHAKARIN AGANFYIVGR DPAGMGHPVE KRDLYDADHG KKVLSMAPGL ERLNILPFRV AAYDKTQGKM AFFDPSRPQD FLFISGTKMR TLAKNNENPP DGFMCPPGGWK VLVDYYESLT PAGNGRLPEV VPV
Source	Mammalian cell
Target Names	APS1
Protein Names	Recommended name: ATP sulfurylase 1, chloroplastic Short name= AtPS1 EC= 2.7.7.4
Expression Region	48-463
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.