



Recombinant Arabidopsis thaliana F-box protein At3g18910 (At3g18910)

Product Code	CSB-EP864773DOA-B
Abbreviation	At3g18910
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9LJ74
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	≥85% (SDS-PAGE)
Sequence	MKTIQEQLPN DLVEEILCRV PATSLRRLRS TCKAWNRLFK GDRILASKHF EKSQKQFRSL SLRNDYRIFP ISFNLHGNP SLELKSELID PHSKNSAAPF EISRVIHCEG LLLCSSQLDE SRVVWNPLT GETRWIRTGD FRQKGRSFDV GYYYQKDKRS WIKSYKLLCY YRGTKYFEIY DFSDSDSWRIL DDIIAPRGS GYSELSVSLK GNTYWFAKGV TEERPRTISL LKFDYFTEKS VPVLLPYQSR RLFQASSLSV VREDKLSVLL QLDQSSKTEI WVTNVIDETT KGAVSWTKVL ALDLSPHLQI GNDGSFFLGE DKKVVMFCEK LIDENKVKDM VYIVGEDNVV TEVGFGVDEM DGCRAVILNY VPSLVQIERA GGNRKRGH
Source	E.coli
Target Names	ETP2
Protein Names	Recommended name: F-box protein At3g18910
Expression Region	1-388
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.