



Recombinant Arabidopsis thaliana Putative F-box protein At3g21130 (At3g21130)

Product Code	CSB-EP864774DOA
Abbreviation	At3g21130
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9LJB9
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MMKKRNTVYL SEDLIVEILS RVSAVSLARL RTTSKRWNAL VKDERLAKKH SAYAPRQSLV ITLIDSRVYL MNVSLQYGIE KVDLSAKLTG QFSLKDPLSN SLEEVDIRNV FHCDGLLCS TKDNRLVWN PCSGETRWIQ PRSSYKVSDI YALGYDNTSS CHKILRMDRS EDRIPIQTEY QVYDFTSKSW LVDGVAGGLF IPSIGTRRRG LSVKGNTYWL ALTEDGPPFD MFLLCFDFST DGFRRLSLPT DTPCTYYDVS LSVTREEQQL CMFKSHGSEL WIATKMESTG AISWSKSHRF HFRIDEMTVL ADQEKPVFVY RFKPKSNMLH IVGEHIYIIQ VDQHSADSKC PYLLTYVPSL VQIQQSI
Source	E.coli
Target Names	At3g21130
Protein Names	Recommended name: Putative F-box protein At3g21130
Expression Region	1-367
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.