



Recombinant Arabidopsis thaliana Putative F-box/kelch-repeat protein At1g15680 (At1g15680)

Product Code	CSB-YP864798DOA
Abbreviation	At1g15680
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9LMR4
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MENREVNWMK SGCKRRIELP EELLAIEVAR LPFISITRFK SVCKGWRS LI ESTYFRHLFV FAHRNSSSSW SLVCGTFGWS VEEMAGFYGC KRYGLPRRLG SYIPPHGLVD KHKIVACTDG LVLLQTVTKR ETFSVGSPVL RQWVQLPPHP WKGISSSVLA IGLVTRVEDG VVMEYKVVCM DIDYRWEVES LILEIYSSLT GTWTRKKVRC PSLIVLSYK RCLSLKMLH WLDTHYRCRS SVGAIVAYDV YADDDEQQFR IIPFPDQKAC FRRAYTTSSG FLVCINKIHL LLRLWRLEEY TSDSGRWQLT QEINLTSFGC DHHYFPVAMH PSETHITYMG NPEKALVSID LKTHKLTLHT KSSAYRDTMV YHYLDTVYVS VEASFDHDFV YTPQFTLPTW MGSVPRSPSI
Source	Yeast
Target Names	At1g15680
Protein Names	Recommended name: Putative F-box/kelch-repeat protein At1g15680
Expression Region	1-410
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.