



Recombinant Arabidopsis thaliana Developmentally regulated G-protein 1 (DRG1)

Product Code	CSB-MP864813DOA
Abbreviation	DRG1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9LQK0
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MGIIERIKEI EAEMARTQKN KATEYHLGQL KAKIAKLRTQ LLEPPKGASG GGEGFEVTKY GHGRVALIGF PSVGKSTLLT MLTGTHSEAA SYEFTTLTCI PGVIHYNDTK IQLLDLPGII EGASEGKGRG RQVIAVAKSS DLVLMVLDAS KSEGHRQILT KELEAVGLRL NKTPPQIYFK KKKTTGGISFN TTAPLTHIDE KLCYQILHEY KIHNAEVLFR ENATVDDFID VIEGNRKYIK CVYVYNKIDV VGIDDVDRLS RQPNSIVISC NLKLNLDRL ARMWDEMGLV RVYSKPQGGQ PDFDEPFVLS SDRGGCTVED FCNHVHRTLK KDMKYALVWG TSTRHNPQNC GLSQHLEDED VVQIVKKKER DEGGRGRFKS HSNAPARIAD REKKAPLKQ
Source	Mammalian cell
Target Names	DRG1
Protein Names	Recommended name: Developmentally regulated G-protein 1 Short name= AtDRG1
Expression Region	1-399
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.