



Recombinant Arabidopsis thaliana Actin-related protein 2 (ARP2)

Product Code	CSB-MP864829DOA
Abbreviation	ARP2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9LSD6
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MDNKNVVVCD NGTGYVKCGF AGENFPTSVF PCVVGRPLL R YEE SLMEQQV KDIVVGETCS ELRHQLDINY PVHNGIVQNW DDMEHVWDHA FYNELKINPS DCKILLTDPPLNPSKNREKM IETMFEKYNF AGVFIQIQAV LTLYAQGLLT GLVIDSGDGV THVVPVVDGY SFPHLTKRMN VAGR HITAYL VDLLSRRGYA MNKTADFETV REIKEKLCYI SYDYKRESQL GLETTILVKN YTLPDGRVIK VGTERFQAPE ALFTPELIDV EGDGMADMVF RCIQEMDIDN RMMMLYQHIVL SGGSTMY PGL PSRLEKEIQD RYLDTVL KGN KDGLK KLR LR IEDPPRRKHM VYLGGA VLAG IMKDAPEFWI NREDYMEEGI NCLNKMSQA
Source	Mammalian cell
Target Names	ARP2
Protein Names	Recommended name: Actin-related protein 2 Alternative name(s): Protein WURM
Expression Region	1-389
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.