



Recombinant Arabidopsis thaliana Calcium-binding protein PBP1 (PBP1)

Product Code	CSB-MP864834DOA
Abbreviation	PBP1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9LSQ6
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MASPKSSTRP NQENQEPQFQ DFFPTMAGKL GGEGLIEEIC KGFELLMDKD KGVITFESLR RNASTVLGLG DLTDDDVR YM INEGDFDRDG ALNQMEFCVL MFRLSPELME ASRCVVTEVI EEEFYRH
Source	Mammalian cell
Target Names	PBP1
Protein Names	Recommended name: Calcium-binding protein PBP1 Alternative name(s): KIC-related protein 2 PINOID-binding protein 1
Expression Region	1-127
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.