



Recombinant Arabidopsis thaliana U-box domain-containing protein 25 (PUB25)

Product Code	CSB-YP864838DOA
Abbreviation	PUB25
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9LT79
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MPRNIEPLDL GIQIPYHFRC PISLELMQDP VTVCTGQTYD RASIESWVSI GNNTTCPVTR APLSDFTLIP NHTLRRLIQE WCVANRSNGV ERIPTPKQPA DPTSVRALLS QASAITGTHV SVRSRAAALR RLRGFARDSD KNRVLIAAHN ATEILIKILF SETTSSSELVS ESLALLVMLP ITEPNQFVSI SSDPGRVEFL TRLLFDSSIE TRVNAAALIE IVSTGTKSAD LKGSISNSES VFEGVLDLLR NPISSRRALK IGIKTLFALC SVKSTRHIAI TAGAPEILID RLAADFDRCD TERALATVEL LCRTPEGCAA FGEHALTVPL LVKTILRVSD RATEYAAGAL LALCTAEERW REEAAGAGVV VQLLLMVQSE CTERAKKKAQ KLLKLLRDSW PDYNSFANS DFGCSSQVVP F
Source	Yeast
Target Names	PUB25
Protein Names	Recommended name: U-box domain-containing protein 25 EC= 6.3.2.- Alternative name(s): Plant U-box protein 25
Expression Region	1-421
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.