



Recombinant Arabidopsis thaliana Probable fructokinase-4 (At3g59480)

Product Code	CSB-MP864920DOA
Abbreviation	At3g59480
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9M1B9
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MAASNGEKSL IVSFGEMLID FVPTVSGVSL ADAPGFIKAP GGAPANVAIA ISRLGGRAAF VGKLGDDDFG HMLAGILKQN GVSAEGINFD TGARTALAFV TLRSDGEREF MFYRNPSADM LLRPDELNLD VIRSAKVFHY GSISLIVEPC RSAHLKAMEV AKEAGALLSY DPNLRLPLWP SKEEAQKQIL SIWDKAEVIK VSDEELMFLT GSDKVDDETA LSLWHSNLKL LLVTLGEKGC RYYTKSFRGS VDPFHVDAVD TTGAGDSFVG ALLCKIVDDR AVLEDEARLR EVLRLANACG AITTTKKGAI PALPTESEVQ SLLKGN
Source	Mammalian cell
Target Names	At3g59480
Protein Names	Recommended name: Probable fructokinase-4 EC= 2.7.1.4
Expression Region	1-326
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.