



Recombinant Pteropus hypomelanus Myc proto-oncogene protein (MYC)

Product Code	CSB-MP865037PLAH
Abbreviation	MYC
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9MZT8
Product Type	Recombinant Protein
Immunogen Species	Pteropus hypomelanus (Island flying fox) (Variable flying fox)
Purity	≥85% (SDS-PAGE)
Sequence	MPLNVSFASR NYDLDYDSVQ PYFYCDEEEN FYHQQQQSEL QPPAPSEDIW KKFELLPTTP LSPSRRLGLC SPSYVAAFAS FSPRDDDDGG GGSFSSADQL EMVTELLGGD MVNQSFICDP DDETFIKNII IQDCMWSGFSA AAKLVSEKL ASYQAARKDG GSRSPARGHS ACSTSSLYLQ DLSAAASECI DPSVVFPYPL NDSSSPKPCA SPDSTAFSPS SDSLLSSAAS SPRASPEPLV LHEETPPTTS SDSEEEQEDE EEIDVVSVEK RQPPAKRSES GSPSAGSHSK PPHSPLVLKR CHVSTHQHNY AAPPSTRKDY PPTKRAKLDG GRVLKQISNN RKCASPRSSD TEENDKRRTH NVLERQRRNE LKRSFFALRD QIPELENNEK APKVVILKKA TAYILAIQAE EQKLISEKDL LRKRREQLKH KLEQLRNSCA
Source	Mammalian cell
Target Names	MYC
Protein Names	Recommended name: Myc proto-oncogene protein Alternative name(s): Proto-oncogene c-Myc Transcription factor p64
Expression Region	1-440
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.