



# Recombinant Pteropus hypomelanus Myc proto-oncogene protein (MYC)

|                          |   |
|--------------------------|---|
| <b>Product Code</b>      | CSB-YP865037PLAH  |
| <b>Abbreviation</b>      | MYC   |
| <b>Storage</b>           | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.  |
| <b>Uniprot No.</b>       | Q9MZT8  |
| <b>Product Type</b>      | Recombinant Protein   |
| <b>Immunogen Species</b> | Pteropus hypomelanus (Island flying fox) (Variable flying fox)  |
| <b>Purity</b>            | ≥85% (SDS-PAGE)   |
| <b>Sequence</b>          | MPLNVSFASR NYDLDYDSVQ PYFYCDEEEN FYHQQQQSEL QPPAPSEDIW<br>KKFELLPTTP LSPSRRLGLC SPSYVAAFAS FSPRDDDDGG GGSFSSADQL<br>EMVTELLGGD MVNQSFICDP DDETFIKNII IQDCMWSGFSA AAKLVSEKL<br>ASYQAARKDG GSRSPARGHS ACSTSSLYLQ DLSAAASECI DPSVVFPYPL<br>NDSSSPKPCA SPDSTAFSPS SDSLLSSAAS SPRASPEPLV LHEETPPTTS<br>SDSEEEQEDE EEIDVVSVEK RQPPAKRSES GSPSAGSHSK PPHSPLVLKR<br>CHVSTHQHNY AAPPSTRKDY PPTKRAKLDS GRVLKQISNN RKCASPRSSD<br>TEENDKRRTH NVLERQRRNE LKRSFFALRD QIPELENNEK APKVVILKKA<br>TAYILAIQAE EQKLISEKDL LRKRREQLKH KLEQLRNSCA |
| <b>Source</b>            | Yeast   |
| <b>Target Names</b>      | MYC   |
| <b>Protein Names</b>     | Recommended name: Myc proto-oncogene protein Alternative name(s): Proto-oncogene c-Myc Transcription factor p64   |
| <b>Expression Region</b> | 1-440   |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.   |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.   |
| <b>Protein Length</b>    | full length protein   |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.   |
| <b>Shelf Life</b>        | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life   |



of lyophilized form is 12 months at -20°C/-80°C.