



Recombinant Human SH2 domain-containing protein 2A (SH2D2A)

Product Code	CSB-MP865088HU
Abbreviation	SH2D2A
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9NP31
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MEFPLAQICP QGSHEAPIPT FSTFQITDMT RRSCQNLGYT AASPQAPEAA SNTGNAERAE EVPGEGSLFL QAETRAWFQK TQAHWLLQHG AAPAWFHGFI TRREAERLLE PKPQGCVLVR FSESAVTFVL TYRSRTCCRH FLLAQLRDGR HVVLGEDSAH ARLQDLLLHY TAHPLSPYGE TLTEPLARQT PEPAGLSLRT EESNFGSKSQ DPNPQYSPII KQGQAPVPMQ KEGAGEKEPS QLLRPKPPIP AKPQLPPEVY TIPVPRHRPA PRPKPSNPIY NEPDEPIAFY AMGRGSPGEA PSNIYVEVED EGLPATLGHP VLRKSWSRPV PGGQNTGGSQ LHESENSVIGQ GPPLPHQPPP AWRHTLPHNL SRQVLQDRGQ AWLPLGPPQ
Source	Mammalian cell
Target Names	SH2D2A
Protein Names	Recommended name: SH2 domain-containing protein 2A Alternative name(s): SH2 domain-containing adapter protein T cell-specific adapter protein Short name= TSA _d VEGF receptor-associated protein
Expression Region	1-389
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.