



Recombinant Human Probable Xaa-Pro aminopeptidase 3 (XPNPEP3)

Product Code	CSB-MP865106HU
Abbreviation	XPNPEP3
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9NQH7
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MPWLLSAPKL VPAVANVRGL SGCMLCSQRR YSLQPVPERR IPNRYLGQPS PFTHPHLLRP GEVTPGLSQV EYALRRHKLM SLIQKEAQQG SGTDQTVVVL SNPTYYSND IPYTFHQDNN FLYLCGFQEP DSILVLQSLP GKQLPSHKAI LFVPRRDPSR ELWDGPRSGT DGAIALTGVD EAYTLEEFQH LLPKMKAEIN MVWYDWMRPS HAQLHSDYMQ PLTEAKAKSK NKVRGVQQLI QRLRLIKSPA EIERMQIAGK LTSQAFIETM FTSKAPVEEA FLYAKFEFEC RARGADILAY PPVVAGGNRS NTLHYVKNNQ LIKDGMVLL DGGCESSCYV SDITRTWPVN GRFTAPQAEL YEAVLEIQRD CLALCFPGTS LENIYSMMLT LIGQKLDLG IMKNIKENNA FKAARKYCPH HVGHYLGMDV HTPDMPRSL PLQPGMVITI EPGIYIPEDD KDAPEKFRGL GVRIEDDVVV TQDSPLILSA DCPKEMNDIE QICSQAS
Source	Mammalian cell
Target Names	XPNPEP3
Protein Names	Recommended name: Probable Xaa-Pro aminopeptidase 3 Short name= X-Pro aminopeptidase 3 EC= 3.4.11.9 Alternative name(s): Aminopeptidase P3 Short name= APP3
Expression Region	1-507
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.