



# Recombinant Human Probable Xaa-Pro aminopeptidase 3 (XPNPEP3)

<b>Product Code</b>	CSB-YP865106HU
<b>Abbreviation</b>	XPNPEP3
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9NQH7
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MPWLLSAPKL VPAVANVRGL SGCMLCSQRR YSLQPVPERR IPNRYLGQPS PFTHPHLLRP GEVTPGLSQV EYALRRHKLM SLIQKEAQQG SGTDQTVVVL SNPTYYSND IPYTFHQDNN FLYLCGFQEP DSILVLQSLP GKQLPSHKAI LFVPRRDPSR ELWDGPRSGT DGAIALTGVD EAYTLEEFQH LLPKMKAEIN MVWYDWMRPS HAQLHSDYMQ PLTEAKAKSK NKVRGVQQLI QRLRLIKSPA EIERMQIAGK LTSQAFIETM FTSKAPVEEA FLYAKFEFEC RARGADILAY PPVVAGGNRS NTLHYVKNNQ LIKDGEVLL DGGCESSCYV SDITRTWPVN GRFTAPQAEL YEAVLEIQRD CLALCFPGTS LENIYSMMLT LIGQKLDLGL IMKNIKENNA FKAARKYCPH HVGHYLGMDV HDTPDMRSL PLQPGMVITI EPGIYIPEDD KDAPEKFRGL GVRIEDDVVV TQDSPLILSA DCPKEMNDIE QICSQAS
<b>Source</b>	Yeast
<b>Target Names</b>	XPNPEP3
<b>Protein Names</b>	Recommended name: Probable Xaa-Pro aminopeptidase 3 Short name= X-Pro aminopeptidase 3 EC= 3.4.11.9 Alternative name(s): Aminopeptidase P3 Short name= APP3
<b>Expression Region</b>	1-507
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.