



Recombinant Human Striatin-4 (STRN4)

Product Code	CSB-MP865120HU
Abbreviation	STRN4
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9NRL3
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MMEERAAAV AAAASSCRPL GSGAGPGPTG AAPVSAPAPG PGPAGKGGGG GGSPGPTAGP EPLSLPGILH FIQHEWARFE AEKARWEAER AELQAQVAFL QGERKGQENL KTDLVRRIM LEYALKQERA KYHKLKFGTD LNQGEKKADV SEQVSNGPVE SVTLENSPLV WKEGRQLLRQ YLEEVGYTDT ILDMRSKRVR SLLGRSLELN GAVEPSEGAP RAPPGPAGLS GGESLLVKQI EEQIKRNAAG KDGKERLGGS VLGQIPFLQN CEDEDSDEDD ELDSVQHKKQ RVKLPSKALV PEMEDEDEED DSEDAINEFD FLGSGEDGEG APDPRRCTVD GSPHELESRR VKLQGILADL RDVDGLPPKV TGPPPGTPQP RPHEDEVFIMD TIGGGEVSLG DLADLTVTND NDLSCDLSDS KDAFKKTWNP KFTLRSHYDG IRSLAFHHSQ SALLTASEDG TLKLWNLQKA VTAKKNAALD VEPIHAFRAH RGPVLAVAMG SNSEYCYSGG ADACIHSWKI PDLSDMPYDG YDPSVLSHVL EGHGDAVWGL AFSPTSQRLA SCSADGTVRI WDPSSSSPAC LCTFPTASEH GVPTSVAFTS TEPAHIVASF RSGDTVLYDM EVGSALLTLE SRGSSGPTQI NQVVSHPNQP LTITAHDDRG IRFLDNRTGK PVHSMVAHLD AVTCLAVDPN GAFLMSGSHD CSLRLWSLDN KTCVQEITAH RKKHEEAIHA VACHPSKALI ASAGADALAK VFV
Source	Mammalian cell
Target Names	STRN4
Protein Names	Recommended name: Striatin-4 Alternative name(s): Zinedin
Expression Region	1-753
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.