



Recombinant Human LanC-like protein 2 (LANCL2)

Product Code	CSB-YP865129HU
Abbreviation	LANCL2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9NS86
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	GETMSKRLK LHLGGAEEME ERAFVNPFDP YEAAAGALLA SGAAEETGCV RPPATTDEPG LPFHQDGKII HNFIRRIQTK IKDLLQQMEE GLKTADPHDC SAYTGWTGIA LLYLQLYRVT CDQTYLLRSL DYVKRTLRLN NGRRVTFLCG DAGPLAVGAV IYHKLRSDCE SQECVTKLLQ LQRSVVCQES DLPDELLYGR AGYLYALLYL NTEIGPGTVC ESAIKEVVNA IIESGKTLR EERKTERCPL LYQWHRKQYV GAAHGMAGIY YMLMQPAKV DQETLTEMVK PSIDYVRHKK FRSGNYPSSL SNETDRLVHW CHGAPGVIHM LMQAYKVFKE EKYLKEAMEC SDVIWQRGLL RKGYGICHT AGNGYSFLSL YRLTQDKKYL YRACKFAEWC LDYGAHGCRIPDRPYSLFEG MAGAIHFLSD VLGPESTRFP AFELDSSKRD
Source	Yeast
Target Names	LANCL2
Protein Names	Recommended name: LanC-like protein 2 Alternative name(s): Testis-specific adriamycin sensitivity protein
Expression Region	2-450
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.