



# Recombinant Human Protein DEPP (DEPP)

<b>Product Code</b>	CSB-MP865135HU
<b>Abbreviation</b>	DEPP
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9NTK1
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MRSRLLLSVA HLPTIRETTE EMLLGGPGQE PPPSPSLDDY VRSISRLAQP TSVLDKATAQ GQPRPPHRPA QACRKGRPAV SLRDITARFS GQQPTLPMAD TVDPLDWLFG ESQEKQPSQR DLPRRTGPSA GLWGPHRQMD SSKPMGAPRG RLCEARMPGH SLARPPQDGQ QSSDLRSWTF GQSAQAMASR HRPRPSSVLR TLYSHLPVIH EL
<b>Source</b>	Mammalian cell
<b>Target Names</b>	DEPP1
<b>Protein Names</b>	Recommended name: Protein DEPP Alternative name(s): Decidual protein induced by progesterone Fasting-induced gene protein Short name= FIG
<b>Expression Region</b>	1-212
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Target Details</b>	The expression of this gene is induced by fasting as well as by progesterone. This protein contains a t-synaptosome-associated protein receptor (SNARE) coiled-coil homology domain and a peroxisomal targeting signal. Production of the encoded protein leads to phosphorylation and activation of the transcription factor ELK1.
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.