



Recombinant Human UPF0609 protein C4orf27 (C4orf27)

Product Code	CSB-EP865162HU-B
Abbreviation	C4orf27
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9NWX4
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MVGGGGKRRP GGEGPQCEKT TDVKKSKFCE ADVSSDLRKE VENHYKLSLP EDFYHFWKFC EELDPEKPSD SLSASLGLQL VGPYDILAGK HKTKKKSTGL NFNLHWRFYD DPPEFQTHI GDNKTQYHMG YFRDSPDEFV VYVGINEAKK NCIIVPNGDN VFAAVKLFLT KKLREITDKK KINLLKNIDE KLTEAARELG YSLEQRTVKM KQRDKKVVTK TFHGAGLVVP VDKNDVGYRE LPETDADLKR ICKTIVEAAS DEERLKAFAP IQEMMTFVQF ANDECDYGMG LELGMDLFCY GSHYFHKVAG QLLPLAYNLL KRNLFAEIII EHLANRSQEN IDQLAA
Source	E.coli
Target Names	HPF1
Protein Names	Recommended name: UPF0609 protein C4orf27
Expression Region	1-346
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.