



Recombinant Human Uncharacterized protein C20orf111 (C20orf111)

Product Code	CSB-MP865164HU
Abbreviation	C20orf111
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9NX31
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MKSEAKDGEE ESLQTAFKKL RVDASGSVAS LSVGEGTGVR APVRTATDDT KPKTTCASKD SWHGSTRKSS RGAVRTQRRR RSKSPVLHPP KFIHCSTIAS SSSSQLKHKS QTDSPDGSSG LGISSPKEFS AGESSTSLDA NHTGAVVEPL RTSVPRLPSE SKKEDSSDAT QVPQASLKAS DLSDFQSVSK LNQGKPCTCI GKECQCKRWH DMEVYSFSG L QSVPLAPER RSTLEDYSQS LHARTLSGSP RSCSEQARVF VDDVTIEDLS GYMEYYLYIP KKMSHMAEMM YT
Source	Mammalian cell
Target Names	OSER1
Protein Names	Recommended name: Uncharacterized protein C20orf111
Expression Region	1-292
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.