



# Recombinant Human Protein AATF (AATF)

<b>Product Code</b>	CSB-MP865175HU
<b>Abbreviation</b>	AATF
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9NY61
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	AGPQPLALQ LEQLLNPRPS EADPEADPEE ATAARVIDRF DEGEDGEGDF LVVGSIRKLA SASLLDTRDKR YCGKTTSRKA WNEDHWEQTL PGSSDEEISD EEGSGDEDSE GLGLEEYDED DLGAAEEQEC GDHRESKKS R SHSAKTPGFS VQSISDFEKF TKGMDLGGSS EEEEEDEESGM EEGDDAEDSQ GESEEDRAGD RNSEDDGVVM TFSSVKVSEE VEKGRAVKNQ IALWDQLLEG RIKLQKALLT TNQLPQPDVF PLFKDKGGPE FSSALKNSHK ALKALLRSLV GLQEELLFQY PDTRYLVDGT KPNAGSEEIS SEDDELVEEK KQRRRRVPAK RKLEMEDYPS FMAKRFADFT VYRNRTLQKW HDKTKLASGK LGKGFGA FER SILTQIDHIL MDKERLLRRT QTKRSVYRVL GKPEPAAQPV PESLPGEPEI LPQAPAN AHL KDLDEEIFDD DDFYHQLLRE LIERKTSSLD PNDQVAMGRQ WLAIQKLR SK IHKKVDRKAS KGRKLRFHVL SKLLSFM API DHTTMNDDAR TELYRSLFGQ LHPPDEGHGD
<b>Source</b>	Mammalian cell
<b>Target Names</b>	AATF
<b>Protein Names</b>	Recommended name: Protein AATF Alternative name(s): Apoptosis-antagonizing transcription factor Rb-binding protein Che-1
<b>Expression Region</b>	2-560
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Target Details</b>	This protein was identified on the basis of its interaction with MAP3K12/DLK, a protein kinase known to be involved in the induction of cell apoptosis. This gene product contains a leucine zipper, which is a characteristic motif of transcription factors, and was shown to exhibit strong transactivation activity when fused to Gal4 DNA binding domain. Overexpression of this gene interfered with MAP3K12 induced apoptosis.
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a



concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

### Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.