



# Recombinant Human Homeobox protein Hox-C10 (HOXC10)

<b>Product Code</b>	CSB-MP865177HU
<b>Abbreviation</b>	HOXC10
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9NYD6
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MTCPRNVTPN SYAEPLAAPG GGERYSRSAG MYMQSGSDFN CGVMRGCGLA PLSKRDEGS SPSLALNTYP SYLSQLDSWG DPKAAYRLEQ PVGRPLSSCS YPPSVKEENV CCMYSAEKRA KSGPEAALYS HPLPESCLGE HEVPVPSYYR ASPSYSALDK TPHCSGANDF EAPFEQRASL NPRAEHLESP QLGGKVSFPE TPKSDSQTPS PNEIKTEQSL AGPKGSPSES EKERAKAADS SPDTSNEAK EEIKAENTTG NWLTAKSGRK KRCPYTKHQT LELEKEFLFN MYLTRERRLE ISKTINLTDR QVKIWFQNRK MKLKKMNREN RIRELTSNFN FT
<b>Source</b>	Mammalian cell
<b>Target Names</b>	HOXC10
<b>Protein Names</b>	Recommended name: Homeobox protein Hox-C10 Alternative name(s): Homeobox protein Hox-3l
<b>Expression Region</b>	1-342
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Target Details</b>	This gene belongs to the homeobox family of genes. The homeobox genes encode a highly conserved family of transcription factors that play an important role in morphogenesis in all multicellular organisms. Mammals possess four similar homeobox gene clusters, HOXA, HOXB, HOXC and HOXD, which are located on different chromosomes and consist of 9 to 11 genes arranged in tandem. This gene is one of several homeobox HOXC genes located in a cluster on chromosome 12. The protein level is controlled during cell differentiation and proliferation, which may indicate this protein has a role in origin activation.
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final



concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

---

**Shelf Life**

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.